

For Research Use Only

ACSF2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84033-1-PBS



Basic Information

Catalog Number:

84033-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9199

GenBank Accession Number:

BC012053

GeneID (NCBI):

80221

UNIPROT ID:

Q96CM8

Full Name:

acyl-CoA synthetase family member

2

Calculated MW:

615 aa, 68 kDa

Purification Method:

Protein A purification

CloneNo.:

241208C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00952-2, ACSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84033-1-PBS. Detection antibody: 84033-2-PBS. Standard: Ag9199. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00952-3, ACSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84033-1-PBS. Detection antibody: 84033-3-PBS. Standard: Ag9199. Range: 2.5-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84033-1-PBS against Human ACSF2 was performed. The affinity constant is below 1 pM.