

For Research Use Only

# UTP18 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84035-4-PBS

## Basic Information

Catalog Number: 84035-4-PBS	GenBank Accession Number: NM_016001	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51096	CloneNo.: 241240G6
Source: Rabbit	UNIPROT ID: Q9Y5J1	
Isotype: IgG	Full Name: UTP18, small subunit (SSU) processome component, homolog (yeast)	
Immunogen Catalog Number: AG34662	Calculated MW: 62 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00941-1, UTP18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84035-1-PBS. Detection antibody: 84035-4-PBS. Standard: Ag34662. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00941-2, UTP18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84035-3-PBS. Detection antibody: 84035-4-PBS. Standard: Ag34662. Range: 0.781-100 ng/mL.