

For Research Use Only

RIPK1-Specific Recombinant antibody, PBS Only (Detector)

Catalog Number: 84036-1-PBS



Basic Information

Catalog Number:

84036-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003804

GeneID (NCBI):

8737

UNIPROT ID:

Q13546

Full Name:

receptor (TNFRSF)-interacting serine-threonine kinase 1

Calculated MW:

76 kDa

Purification Method:

Protein A purification

CloneNo.:

241099C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00944-1, RIPK1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84036-2-PBS. Detection antibody: 84036-1-PBS. Standard: Ag31943. Range: 1.563-100 ng/mL.