

For Research Use Only

GPAM Recombinant antibody, PBS Only (Detector)

Catalog Number: 84054-1-PBS



Basic Information

Catalog Number:

84054-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3120

GenBank Accession Number:

BC030783

GeneID (NCBI):

57678

UNIPROT ID:

Q9HCL2

Full Name:

glycerol-3-phosphate acyltransferase,
mitochondrial

Calculated MW:

828 aa, 94 kDa

Purification Method:

Protein A purification

CloneNo.:

241144B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00994-2, GPAM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84054-5-PBS. Detection antibody: 84054-1-PBS. Standard: Ag3120. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84054-1-PBS against Human GPAM was performed. The affinity constant is 0.846 nM.