

For Research Use Only

MRP3/ABCC3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84057-1-PBS



Basic Information

Catalog Number:

84057-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19786

GenBank Accession Number:

BC137348

GeneID (NCBI):

8714

UNIPROT ID:

O15438

Full Name:

ATP-binding cassette, sub-family C
(CFTR/MRP), member 3

Calculated MW:

1527 aa, 169 kDa

Purification Method:

Protein A purification

CloneNo.:

241209B9

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00973-2, Human MRP3/ABCC3 Recombinant Matched Antibody Pair - PBS only. 84057-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19786. 84057-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



The mean MRP3/ABCC3 concentration was determined to be 1.31 ng/mL in HeLa cell extract based on a 2.0 mg/mL extract load.



Cytometric bead array standard curve of MP00973-1, MRP3/ABCC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84057-1-PBS. Detection antibody: 84057-2-PBS. Standard: Ag19786. Range: 0.313-40 ng/mL