

For Research Use Only

MRP3/ABCC3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84057-2-PBS



Basic Information

Catalog Number:

84057-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19786

GenBank Accession Number:

BC137348

GeneID (NCBI):

8714

UNIPROT ID:

O15438

Full Name:

ATP-binding cassette, sub-family C
(CFTR/MRP), member 3

Calculated MW:

1527 aa, 169 kDa

Purification Method:

Protein A purification

CloneNo.:

241209E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00973-1, MRP3/ABCC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84057-1-PBS. Detection antibody: 84057-2-PBS. Standard: Ag19786. Range: 0.313-40 ng/mL.