

For Research Use Only

SRP19 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84061-1-PBS



Basic Information

Catalog Number:

84061-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8903

GenBank Accession Number:

BC010947

GeneID (NCBI):

6728

UNIPROT ID:

P09132

Full Name:

signal recognition particle 19kDa

Calculated MW:

144 aa, 16 kDa

Purification Method:

Protein A purification

CloneNo.:

241219C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00977-1, SRP19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84061-1-PBS. Detection antibody: 84061-3-PBS. Standard: Ag8903. Range: 0.195-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84061-1-PBS against Human SRP19 was performed. The affinity constant is below 1 pM.