

For Research Use Only

# SRP19 Recombinant antibody, PBS Only

Catalog Number: 84061-2-PBS



## Basic Information

Catalog Number:

84061-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8903

GenBank Accession Number:

BC010947

GeneID (NCBI):

6728

UNIPROT ID:

P09132

Full Name:

signal recognition particle 19kDa

Calculated MW:

144 aa, 16 kDa

Purification Method:

Protein A purification

CloneNo.:

241219D12

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

