

For Research Use Only

LRRK2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84062-1-PBS



Basic Information

Catalog Number:	84062-1-PBS	GenBank Accession Number:	NM_198578	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	120892	CloneNo.:	241258D3
Source:	Rabbit	UNIPROT ID:	Q5S007		
Isotype:	IgG	Full Name:	leucine-rich repeat kinase 2		
Immunogen Catalog Number:	AG33713	Calculated MW:	286 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00980-1, LRRK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84062-1-PBS. Detection antibody: 84062-2-PBS. Standard: Ag33713. Range: 0.625-80 ng/mL