

For Research Use Only

# LRRK2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84062-2-PBS



## Basic Information

Catalog Number: 84062-2-PBS	GenBank Accession Number: NM_198578	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 120892	CloneNo.: 241258D8
Source: Rabbit	UNIPROT ID: Q5S007	
Isotype: IgG	Full Name: leucine-rich repeat kinase 2	
Immunogen Catalog Number: AG33713	Calculated MW: 286 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00980-1, LRRK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84062-1-PBS. Detection antibody: 84062-2-PBS. Standard: Ag33713. Range: 0.625-80 ng/mL