

For Research Use Only

# TICAM2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84064-2-PBS



## Basic Information

Catalog Number:

84064-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_021649

GeneID (NCBI):

353376

UNIPROT ID:

Q86XR7

Full Name:

toll-like receptor adaptor molecule 2

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

241186D8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00979-1, TICAM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84064-2-PBS. Detection antibody: 84064-4-PBS. Standard: Ag17325. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00979-3, TICAM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84064-2-PBS. Detection antibody: 84064-1-PBS. Standard: Ag17325. Range: 0.625-80 ng/mL.