

For Research Use Only

Clusterin Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84067-2-PBS



Basic Information

Catalog Number:

84067-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC010514

GeneID (NCBI):

1191

UNIPROT ID:

P10909

Full Name:

clusterin

Calculated MW:

449 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

241277E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00986-1, Clusterin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84067-2-PBS. Detection antibody: 84067-3-PBS. Standard: Eg0225. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00986-2, Clusterin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84067-1-PBS. Detection antibody: 84067-2-PBS. Standard: Eg0225. Range: 0.625-40 ng/mL.