

For Research Use Only

Leptin Recombinant antibody, PBS Only (Detector)

Catalog Number: 84072-3-PBS



Basic Information

Catalog Number:

84072-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC060830

GeneID (NCBI):

3952

UNIPROT ID:

P41159

Full Name:

leptin

Calculated MW:

167 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

240720D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00993-3, LEP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84072-2-PBS. Detection antibody: 84072-3-PBS. Standard: Eg0834. Range: 1.25-20 ng/mL.