

For Research Use Only

# FKBP2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84077-2-PBS

## Basic Information

Catalog Number: 84077-2-PBS	GenBank Accession Number: BC003384	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2286	CloneNo.: 241166C10
Source: Rabbit	UNIPROT ID: P26885	
Isotype: IgG	Full Name: FK506 binding protein 2, 13kDa	
Immunogen Catalog Number: AG2255	Calculated MW: 142 aa, 16 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00984-1, FKBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84077-1-PBS. Detection antibody: 84077-2-PBS. Standard: Ag2255. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00984-2, FKBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84077-5-PBS. Detection antibody: 84077-2-PBS. Standard: Ag2255. Range: 0.781-100 ng/mL.