

For Research Use Only

SEMA3B Recombinant antibody, PBS Only (Capture)

Catalog Number: 84080-1-PBS



Basic Information

Catalog Number:

84080-1-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25163

GenBank Accession Number:

BC009113

GeneID (NCBI):

7869

UNIPROT ID:

Q13214

Full Name:

sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3B

Purification Method:

Protein A purification

CloneNo.:

241171B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00995-1, SEMA3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84080-1-PBS. Detection antibody: 84080-2-PBS. Standard: Ag25163. Range: 0.625-80 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84080-1-PBS against Human SEMA3B was performed. The affinity constant is below 1 pM.