

For Research Use Only

TRP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84082-1-PBS



Basic Information

Catalog Number:

84082-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10166

GenBank Accession Number:

BC052608

GeneID (NCBI):

7306

UNIPROT ID:

P17643

Full Name:

tyrosinase-related protein 1

Calculated MW:

537 aa, 61 kDa

Purification Method:

Protein A purification

CloneNo.:

241242B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00997-1, TRP1/TYRP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84082-1-PBS. Detection antibody: 84082-2-PBS. Standard: Ag10166. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00997-2, TRP1/TYRP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84082-1-PBS. Detection antibody: 84082-3-PBS. Standard: Ag10166. Range: 0.625-80 ng/mL.