

For Research Use Only

CD209 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84095-2-PBS



Basic Information

Catalog Number:

84095-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0062

GenBank Accession Number:

BC110615

GeneID (NCBI):

30835

UNIPROT ID:

Q9NNX6

Full Name:

CD209 molecule

Calculated MW:

404 aa, 45.7 kDa

Purification Method:

Protein A purification

CloneNo.:

240907A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01013-1, CD209 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84095-1-PBS. Detection antibody: 84095-2-PBS. Standard: Eg0062. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01013-2, CD209 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84095-3-PBS. Detection antibody: 84095-2-PBS. Standard: Eg0062. Range: 0.313-20 ng/mL