

For Research Use Only

RNASEH2B Recombinant antibody, PBS Only (Capture)



Catalog Number: 84097-1-PBS

Basic Information

Catalog Number: 84097-1-PBS	GenBank Accession Number: BC036744	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79621	CloneNo.: 241131C10
Source: Rabbit	UNIPROT ID: Q5TBB1	
Isotype: IgG	Full Name: ribonuclease H2, subunit B	
Immunogen Catalog Number: AG33669	Calculated MW: 331 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01015-1, RNASEH2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84097-1-PBS. Detection antibody: 84097-5-PBS. Standard: Ag33669. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01015-3, RNASEH2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84097-1-PBS. Detection antibody: 84097-3-PBS. Standard: Ag33669. Range: 0.781-50 ng/mL