

For Research Use Only

ACMSD Recombinant antibody, PBS Only (Capture)

Catalog Number: 84100-2-PBS



Basic Information

Catalog Number:

84100-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25212

GenBank Accession Number:

BC016018

GeneID (NCBI):

130013

UNIPROT ID:

Q8TDX5

Full Name:

aminocarboxymuconate
semialdehyde decarboxylase

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

241077F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01022-1, ACMSD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84100-2-PBS. Detection antibody: 84100-1-PBS. Standard: Ag25212. Range: 2.5-80 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84100-2-PBS against Human ACMSD was performed. The affinity constant is 0.886 nM.