

For Research Use Only

TSEN34 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84102-1-PBS

Basic Information

Catalog Number: 84102-1-PBS	GenBank Accession Number: BC020805	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79042	CloneNo.: 241155B7
Source: Rabbit	UNIPROT ID: Q9BSV6	
Isotype: IgG	Full Name: tRNA splicing endonuclease 34 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG35120	Calculated MW: 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01036-2, TSEN34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84102-1-PBS. Detection antibody: 84102-4-PBS. Standard: Ag35120. Range: 0.313-40 ng/mL