

For Research Use Only

TSEN34 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 84102-2-PBS

Basic Information

Catalog Number:

84102-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35120

GenBank Accession Number:

BC020805

GeneID (NCBI):

79042

UNIPROT ID:

Q9BSV6

Full Name:

tRNA splicing endonuclease 34 homolog (S. cerevisiae)

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

241155D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01036-1, TSEN34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84102-3-PBS. Detection antibody: 84102-2-PBS. Standard: Ag35120. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01036-3, TSEN34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84102-3-PBS. Detection antibody: 84102-4-PBS. Standard: Ag35120. Range: 0.625-40 ng/mL