

For Research Use Only

# CHRNA2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84116-4-PBS



## Basic Information

Catalog Number:

84116-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12208

GenBank Accession Number:

BC075041

GeneID (NCBI):

1141

UNIPROT ID:

P17787

Full Name:

cholinergic receptor, nicotinic, beta 2 (neuronal)

Calculated MW:

502 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

241374G9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01019-2, CHRNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84116-1-PBS. Detection antibody: 84116-4-PBS. Standard: Ag12208. Range: 0.625-80 ng/mL.