

For Research Use Only

# PDLIM4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84117-4-PBS

## Basic Information

Catalog Number: 84117-4-PBS	GenBank Accession Number: BC016765	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8572	CloneNo.: 241401H7
Source: Rabbit	UNIPROT ID: P50479	
Isotype: IgG	Full Name: PDZ and LIM domain 4	
Immunogen Catalog Number: AG13424	Calculated MW: 35 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01021-1, PDLIM4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84117-4-PBS. Detection antibody: 84117-1-PBS. Standard: Ag13424. Range: 0.781-100 ng/mL.