

For Research Use Only

SCF Recombinant antibody, PBS Only (Detector)



Catalog Number: 84118-1-PBS

Basic Information

Catalog Number:

84118-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074725

GeneID (NCBI):

4254

UNIPROT ID:

P21583

Full Name:

KIT ligand

Calculated MW:

273 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

241365D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01028-2, SCF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84118-2-PBS. Detection antibody: 84118-1-PBS. Standard: Eg0112. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP01028-3, SCF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84118-4-PBS. Detection antibody: 84118-1-PBS. Standard: Eg0112. Range: 0.625-40 ng/mL