

For Research Use Only

SCF Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 84118-9-PBS



Basic Information

Catalog Number:

84118-9-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1024

GenBank Accession Number:

GeneID (NCBI):

4254

UNIPROT ID:

P21583

Full Name:

KIT ligand

Purification Method:

Protein A purification

CloneNo.:

250446B8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of seven individual healthy human donors was measured. The human SCF concentration of detected samples was determined to be 4.91 ng/mL with a range of 2.40 - 14.68 ng/mL.



Sandwich ELISA standard curve of MP01028-4, Human SCF Recombinant Matched Antibody Pair - PBS only. 84118-9-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1024. 84118-8-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL