

For Research Use Only

# ZNF350 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84121-3-PBS

## Basic Information

Catalog Number:

84121-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23982

GenBank Accession Number:

BC009921

GeneID (NCBI):

59348

UNIPROT ID:

Q9GZX5

Full Name:

zinc finger protein 350

Purification Method:

Protein A purification

CloneNo.:

241225B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01034-2, ZNF350 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84121-3-PBS. Detection antibody: 84121-1-PBS. Standard: Ag23982. Range: 0.781-100 ng/mL.