

For Research Use Only

# SYT9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84123-1-PBS



## Basic Information

Catalog Number:

84123-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3520

GenBank Accession Number:

BC029605

GeneID (NCBI):

143425

UNIPROT ID:

Q86556

Full Name:

synaptotagmin IX

Calculated MW:

491 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

241407D1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01023-1, SYT9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84123-2-PBS. Detection antibody: 84123-1-PBS. Standard: Ag3520. Range: 0.781-100 ng/mL.