

For Research Use Only

# HESX1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84125-4-PBS



## Basic Information

Catalog Number:

84125-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12348

GenBank Accession Number:

BC093979

GeneID (NCBI):

8820

UNIPROT ID:

Q9UBX0

Full Name:

HESX homeobox 1

Calculated MW:

185 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

241430G8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01027-1, HESX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84125-4-PBS. Detection antibody: 84125-3-PBS. Standard: Ag12348. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01027-3, HESX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84125-2-PBS. Detection antibody: 84125-4-PBS. Standard: Ag12348. Range: 0.781-100 ng/mL.