

For Research Use Only

# HTR2A Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84126-1-PBS



## Basic Information

Catalog Number:

84126-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24056

GenBank Accession Number:

BC074848

GeneID (NCBI):

3356

UNIPROT ID:

P28223

Full Name:

5-hydroxytryptamine (serotonin)  
receptor 2A

Calculated MW:

471 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

241440C12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01029-2, 5HT2A Receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84126-1-PBS. Detection antibody: 84126-3-PBS. Standard: Ag24056. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01029-3, 5HT2A Receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84126-2-PBS. Detection antibody: 84126-1-PBS. Standard: Ag24056. Range: 0.625-40 ng/mL.