

For Research Use Only

# RAB9A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84129-1-PBS



## Basic Information

Catalog Number:

84129-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1969

GenBank Accession Number:

BC017265

GeneID (NCBI):

9367

UNIPROT ID:

P51151

Full Name:

RAB9A, member RAS oncogene family

Calculated MW:

201 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

241372A6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01043-3, RAB9A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84129-4-PBS. Detection antibody: 84129-1-PBS. Standard: Ag1969. Range: 0.313-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84129-1-PBS against Human RAB9A was performed. The affinity constant is 0.453 nM.