

For Research Use Only

IMUP Recombinant antibody, PBS Only (Capture)

Catalog Number: 84148-4-PBS



Basic Information

Catalog Number:	84148-4-PBS	GenBank Accession Number:	BC060319	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	64073	CloneNo.:	241220D10
Source:	Rabbit	UNIPROT ID:	Q9GZP8		
Isotype:	IgG	Full Name:	chromosome 19 open reading frame		
Immunogen Catalog Number:	AG22218	33			
		Calculated MW:	106 aa, 11 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01103-2, C19orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84148-4-PBS. Detection antibody: 84148-3-PBS. Standard: Ag22218. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01103-3, C19orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84148-4-PBS. Detection antibody: 84148-1-PBS. Standard: Ag22218. Range: 0.313-40 ng/mL.