

For Research Use Only

# ARHGEF11 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84149-2-PBS



## Basic Information

Catalog Number:

84149-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33104

GenBank Accession Number:

GeneID (NCBI):

9826

UNIPROT ID:

O15085

Full Name:

Rho guanine nucleotide exchange factor (GEF) 11

Calculated MW:

168 kDa

Purification Method:

Protein A purification

CloneNo.:

241427A12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Cytometric bead array standard curve of MP01063-1, ARHGEF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84149-2-PBS. Detection antibody: 84149-1-PBS. Standard: Ag33104. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01063-3, ARHGEF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84149-4-PBS. Detection antibody: 84149-2-PBS. Standard: Ag33104. Range: 0.313-40 ng/mL.



Biolayer Interferometry (BLI) kinetic assay of 84149-2-PBS against Human ARHGEF11 was performed. The affinity constant is 0.262 nM.