

For Research Use Only

FAM136A Recombinant antibody, PBS Only (Capture)

Catalog Number: 84150-3-PBS



Basic Information

Catalog Number:

84150-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34058

GenBank Accession Number:

BC014975

GeneID (NCBI):

84908

UNIPROT ID:

Q96C01

Full Name:

family with sequence similarity 136,
member A

Purification Method:

Protein A purification

CloneNo.:

241433F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01066-1, FAM136A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84150-3-PBS. Detection antibody: 84150-4-PBS. Standard: Ag34058. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01066-2, FAM136A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84150-3-PBS. Detection antibody: 84150-1-PBS. Standard: Ag34058. Range: 0.313-40 ng/mL