

For Research Use Only

# FAM136A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84150-4-PBS



## Basic Information

Catalog Number:

84150-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34058

GenBank Accession Number:

BC014975

GeneID (NCBI):

84908

UNIPROT ID:

Q96C01

Full Name:

family with sequence similarity 136, member A

Purification Method:

Protein A purification

CloneNo.:

241433F10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01066-1, FAM136A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84150-3-PBS. Detection antibody: 84150-4-PBS. Standard: Ag34058. Range: 0.313-40 ng/mL.