

For Research Use Only

# RENBP Recombinant antibody, PBS Only (Capture)

Catalog Number: 84152-2-PBS



## Basic Information

Catalog Number:

84152-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35056

GenBank Accession Number:

BC015558

GeneID (NCBI):

5973

UNIPROT ID:

P51606

Full Name:

renin binding protein

Calculated MW:

49 kDa

Purification Method:

Protein A purification

CloneNo.:

241142E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP01073-1, RENBP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84152-2-PBS. Detection antibody: 84152-1-PBS. Standard: Ag35056. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01073-2, RENBP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84152-2-PBS. Detection antibody: 84152-3-PBS. Standard: Ag35056. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01073-3, RENBP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84152-2-PBS. Detection antibody: 84152-4-PBS. Standard: Ag35056. Range: 0.781-100 ng/mL.