

For Research Use Only

# AMH Recombinant antibody, PBS Only (Detector)

Catalog Number: 84153-1-PBS



## Basic Information

Catalog Number:

84153-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5670

GenBank Accession Number:

BC049194

GeneID (NCBI):

268

UNIPROT ID:

P03971

Full Name:

anti-Mullerian hormone

Calculated MW:

59 kDa

Purification Method:

Protein A purification

CloneNo.:

241148F7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP01075-1, AMH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84153-2-PBS. Detection antibody: 84153-1-PBS. Standard: Ag5670. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01075-3, AMH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84153-4-PBS. Detection antibody: 84153-1-PBS. Standard: Ag5670. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84153-1-PBS against Human AMH was performed. The affinity constant is 4.59 nM.