

For Research Use Only

PPP2CA Recombinant antibody, PBS Only (Detector)

Catalog Number: 84155-3-PBS



Basic Information

Catalog Number:

84155-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4294

GenBank Accession Number:

BC031696

GeneID (NCBI):

5515

UNIPROT ID:

P67775

Full Name:

protein phosphatase 2 (formerly 2A), catalytic subunit, alpha isoform

Calculated MW:

309 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

241112G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01058-1, PPP2CA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84155-1-PBS. Detection antibody: 84155-3-PBS. Standard: Ag4294. Range: 0.781-100 ng/mL.