

For Research Use Only

# ABCC10 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84164-2-PBS



## Basic Information

Catalog Number:

84164-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33719

GenBank Accession Number:

NM\_001198934.1

GeneID (NCBI):

89845

UNIPROT ID:

Q5T3U5

Full Name:

ATP-binding cassette, sub-family C (CFTR/MRP), member 10

Calculated MW:

161kd

Purification Method:

Protein A purification

CloneNo.:

241429E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01067-2, ABCC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84164-3-PBS. Detection antibody: 84164-2-PBS. Standard: Ag33719. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84164-2-PBS against Human ABCC10 was performed. The affinity constant is below 1 pM.