

For Research Use Only

Mouse CD6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84166-4-PBS



Basic Information

Catalog Number:

84166-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009852.3

GeneID (NCBI):

12511

UNIPROT ID:

Q61003-1

Full Name:

CD6 antigen

Calculated MW:

72kDa

Purification Method:

Protein A purification

CloneNo.:

241471A7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01099-3, Mouse CD6 Recombinant Matched Antibody Pair - PBS only. 84166-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1401. 84166-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



Serum of twelve mice was measured. The mouse CD6 concentration of detected samples was determined to be 2.83 ng/mL with a range of 1.10 - 4.98 ng/mL.