

For Research Use Only

TRMT12 Recombinant antibody

Catalog Number: 84167-3-RR

Featured Product



Basic Information

Catalog Number:

84167-3-RR

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9258

GenBank Accession Number:

BC011713

GeneID (NCBI):

55039

UNIPROT ID:

Q53H54

Full Name:

tRNA methyltransferase 12 homolog
(S. cerevisiae)

Calculated MW:

448 aa, 50 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein A purification

CloneNo.:

241483H1

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, MCF-7 cells, K-562 cells, U2OS cells, U-521 cells, MDA-MB-231 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84167-3-RR (TRMT12 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



WB result of TRMT12 antibody (84167-3-RR; 1:1200; incubated at room temperature for 1.5 hours) with sh-Control and sh-TRMT12 transfected MCF-7 cells.