

For Research Use Only

# CHRNA4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84168-2-PBS



## Basic Information

Catalog Number:

84168-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16584

GenBank Accession Number:

BC096290

GeneID (NCBI):

1137

UNIPROT ID:

P43681

Full Name:

cholinergic receptor, nicotinic, alpha 4

Calculated MW:

627 aa, 70 kDa

Purification Method:

Protein A purification

CloneNo.:

241492B12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01101-2, CHRNA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84168-2-PBS. Detection antibody: 84168-4-PBS. Standard: Ag16584. Range: 0.781-100 ng/mL