

For Research Use Only

PRKAG2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84169-1-PBS



Basic Information

Catalog Number:

84169-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3214

GenBank Accession Number:

BC020540

GeneID (NCBI):

51422

UNIPROT ID:

Q9UGJ0

Full Name:

protein kinase, AMP-activated,
gamma 2 non-catalytic subunit

Calculated MW:

569 aa, 63 kDa

Purification Method:

Protein A purification

CloneNo.:

241253A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01102-1, PRKAG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84169-1-PBS. Detection antibody: 84169-2-PBS. Standard: Ag3214. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84169-1-PBS against Human PRKAG2 was performed. The affinity constant is below 1 pM.