

For Research Use Only

# GEMIN5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84171-1-PBS



## Basic Information

Catalog Number:	84171-1-PBS	GenBank Accession Number:	BC113614	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	25929	CloneNo.:	241205A2
Source:	Rabbit	UNIPROT ID:	Q8TEQ6		
Isotype:	IgG	Full Name:	gem (nuclear organelle) associated protein 5		
Immunogen Catalog Number:	AG19220	Calculated MW:	1508 aa, 169 kDa		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01070-1, GEMIN5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84171-1-PBS. Detection antibody: 84171-3-PBS. Standard: Ag19220. Range: 0.781-100 ng/mL