

For Research Use Only

GEMIN5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84171-4-PBS



Basic Information

Catalog Number:

84171-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19220

GenBank Accession Number:

BC113614

GeneID (NCBI):

25929

UNIPROT ID:

Q8TEQ6

Full Name:

gem (nuclear organelle) associated protein 5

Calculated MW:

1508 aa, 169 kDa

Purification Method:

Protein A purification

CloneNo.:

241205E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01070-2, GEMIN5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84171-4-PBS. Detection antibody: 84171-2-PBS. Standard: Ag19220. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01070-3, GEMIN5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84171-4-PBS. Detection antibody: 84171-3-PBS. Standard: Ag19220. Range: 0.781-100 ng/mL