

For Research Use Only

# PCSK9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84172-3-PBS



## Basic Information

Catalog Number:

84172-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_174936.4

GeneID (NCBI):

255738

UNIPROT ID:

Q8NBP7

Full Name:

proprotein convertase  
subtilisin/kexin type 9

Calculated MW:

74 kDa

Purification Method:

Protein A purification

CloneNo.:

241395E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01071-3, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-4-PBS. Detection antibody: 84172-3-PBS. Standard: Eg0362. Range: 1.563-100 ng/mL.