

For Research Use Only

# DPP3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84173-2-PBS



## Basic Information

Catalog Number:

84173-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1022

GenBank Accession Number:

BC007221

GeneID (NCBI):

10072

UNIPROT ID:

Q9NY33

Full Name:

dipeptidyl-peptidase 3

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

241360D2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01074-1, DPP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84173-1-PBS. Detection antibody: 84173-2-PBS. Standard: Ag1022. Range: 1.563-200 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84173-2-PBS against Human DPP3 was performed. The affinity constant is below 1 pM.