

For Research Use Only

# ABLIM2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84177-1-PBS



## Basic Information

Catalog Number:

84177-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18198

GenBank Accession Number:

BC122567

GeneID (NCBI):

84448

UNIPROT ID:

Q6H8Q1

Full Name:

actin binding LIM protein family,  
member 2

Calculated MW:

611 aa, 68 kDa

Purification Method:

Protein A purification

CloneNo.:

241450B12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01092-1, ABLIM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84177-2-PBS. Detection antibody: 84177-1-PBS. Standard: Ag18198. Range: 0.313-40 ng/mL