

For Research Use Only

EPS15 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84179-3-PBS



Basic Information

Catalog Number: 84179-3-PBS	GenBank Accession Number: NM_001981.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2060	CloneNo.: 241485G5
Source: Rabbit	UNIPROT ID: P42566	
Isotype: IgG	Full Name: epidermal growth factor receptor pathway substrate 15	
Immunogen Catalog Number: AG35282	Calculated MW: 99 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01095-3, EPS15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84179-4-PBS. Detection antibody: 84179-3-PBS. Standard: Ag35282. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84179-3-PBS against Human EPS15 was performed. The affinity constant is 4.97 nM.