

For Research Use Only

HAVCR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84183-1-PBS



Basic Information

Catalog Number:

84183-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC020843

GeneID (NCBI):

84868

UNIPROT ID:

Q8TDQ0

Full Name:

hepatitis A virus cellular receptor 2

Calculated MW:

301 aa, 33 kDa

Purification Method:

Protein A purification

CloneNo.:

241421A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01083-1, TIM3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84183-4-PBS. Detection antibody: 84183-1-PBS. Standard: RP02340. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01083-3, TIM3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84183-2-PBS. Detection antibody: 84183-1-PBS. Standard: RP02340. Range: 0.781-100 ng/mL