

For Research Use Only

PIAS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84189-2-PBS



Basic Information

Catalog Number:	84189-2-PBS	GenBank Accession Number:	BC015190	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	9063	CloneNo.:	241474C11
Source:	Rabbit	UNIPROT ID:	O75928		
Isotype:	IgG	Full Name:	protein inhibitor of activated STAT, 2		
Immunogen Catalog Number:	AG9032	Calculated MW:	63 kDa, 69 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01062-2, PIAS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84189-4-PBS. Detection antibody: 84189-2-PBS. Standard: Ag9032. Range: 0.625-80 ng/mL